



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING


PRODUCT NO.: 6324



Mixer, Mazerustar, MD Mixing Adapter

DESCRIPTION: This adapter allows mixing directly inside the Medisca MD Pumps (30 mL, 50 mL, 80 mL and 100 mL) using the Mazerustar KK-300 and the KK-400.

ADDITIONAL ACCESSORIES:

-  6129: Mazerustar KK-400W Spacer for Mixing Adapters (**REQUIRED FOR KK-400 and KK-1000**)

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM