



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6326



ELECTRONIC MORTAR AND PESTLE JAR MIXING ADAPTER (100 mL)

For Mazerustar KK-300SS and KK-400

DESCRIPTION: To be used with the Unguator 100mL Electronic Mortar and Pestle jars. It will hold the 100mL Electronic Mortar and Pestle jar in 6363-01 (Mixer, Mazerustar KK-300SS, Standard Container for Electronic Mortar and Pestle Adapter Insert) to mix in the Mazerustar.