



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6327



MIXING ADAPTER (15 – 50 mL)

For Mazerustar KK-300SS and KK-400

DESCRIPTION: To be used with 15 mL-50 mL Electronic Mortar and Pestle jars. It will hold 15 mL-50 mL Electronic Mortar and Pestle jars in 6363-01 (Mixer, Mazerustar KK-300SS, Standard Container for Electronic Mortar and Pestle Adapter Insert) to mix in the Mazerustar. Must use with 6363-01.