



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6348



## DUAL GRADUATED BOTTLE (3 oz)

Amber, Polypropylene, w/Child Resistant Cap

**DESCRIPTION:** This bottle is ideal for liquids, solids and powders and is clearly marked with both ounces and milliliters which helps to eliminate waste by providing accurate measurements directly on the container.

### TECHNICAL SPECIFICATIONS:

Capacity (oz)	3	
Neck size	33-400	
Material	Polypropylene	
Color	Bottle	Amber
	Cap	White