




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6402



WASTE CHUTE STOPPER, LABCONCO

For XPert and RXPert Enclosures

DESCRIPTION: This waste chute stopper helps to prevent contaminants from escaping from the enclosure through the side wall, when the waste chute is not in use.

TECHNICAL SPECIFICATIONS:

Diameter (in)	5
Shipping Weight (kg / lb)	0.5 / 1
Material	Solid Rubber