



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6408**



### **JAR and CAP (8 oz, 70–400)**

Polypropylene w/Unlined and Ribbed Polypropylene Cap, White

**DESCRIPTION:** This jar can be used for various purposes as required.

#### **TECHNICAL SPECIFICATIONS:**

Dimension, Ø x H (mm)	64 x 86
Capacity / Volume (mL / oz)	240 / 8
Style	Jar
Material	Polypropylene
Neck Finish	70-400
Attached Cap	White Polypropylene Ribbed Unlined Cap