



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6416**



**STAINLESS STEEL LINER**

316 Stainless Steel, For Mazerustar KK-300SS & KK-400W

**DESCRIPTION:** This stainless steel liner goes in the standard mixing container when using the Mazerustar for micronizing purposes.

**TECHNICAL SPECIFICATIONS:**

Material	316 Stainless Steel
Maximum Capacity (mL)	250
Maximum Suggested Capacity (mL)	130