



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6451



LED LIGHT KIT, 4' (Flow Sciences)

For 48" × 30" Enclosures

DESCRIPTION: This LED light helps to keep the work area in the enclosure illuminated. It has sophisticated touch controls and functions such as intensity memory which allows the user to get a preset level of lighting, every time. The angle lighting is adjustable and provides maximum work space visibility

TECHNICAL SPECIFICATIONS:

Lifetime	About 50000 hours
Input Power	575 ma at 12 Volts DC
Lumens	Max of 675 lumens
Beam Angle	120 Degree Beam Angle