



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6539**



## **REPLACEMENT HEPA FILTER, FAN SIZE: 18" × 36" (Flow Sciences)**

**DESCRIPTION:** The HEPA filter protects the fan from particle contamination. It is 4-inches high for extended life. Using the proper change method, Personnel Protective Equipment (PPE) and properly trained personnel, the filters can be changed flawlessly, with either the Bag-out procedure, with no cross-contamination to either the person or the surrounding.