




**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6670**



**PRE-FILTER, AIR SCIENCE<sup>®</sup> (Pack of 6)**  
For Purair<sup>®</sup> Basic Ductless Fume Hood

**DESCRIPTION:** The Purair<sup>®</sup> Basic Ductless Fume Hood Pre-Filter is an electrostatic pre-filter that is 99.5% effective. It is accessible from inside the chamber to contain the release of any particulates that it traps. The Pre-Filter can be changed while the unit is operating to prevent operator exposure to chemical vapors.

Use alongside the main filter to remove contaminants from the air.

**TECHNICAL SPECIFICATION:**

Nominal weight (kg / lb)	0.45 / 1.00
--------------------------	-------------