



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING






PRODUCT NO.: 6683



Clave CSTD Bag Spike Model CS-10

DESCRIPTION: This drug transfer device, which can be used on any solution container, enables safe injection of liquids into IV bags for safe administration. It contains a passive self-sealing mechanism that cannot be deactivated by the user and remains protective from preparation through to its disposal.

FEATURES:

-  No natural rubber latex
-  Non-pyrogenic fluid path
-  Classified as a CSTD (Closed System Transfer Device) under NIOSH guidelines
-  Needle-free system assures compliance with safe handling policies
-  Easy-to-use system requires no assembly of components