



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6729**



**TELESCOPING BASE STAND w/CASTERS, LABCONCO (6 ft)  
W 72.0" x D 29.0"**

**DESCRIPTION:** The telescoping base stand with casters may be adjusted to eight height positions prior to installation. Telescoping legs are adjustable in 1-inch increments from 27.5 – 34 inches. Casters provide mobility for short distances, such as to gain access to the area behind the stand for service or maintenance. Toe-locking, non-marking polyurethane casters, 5-inch diameter, have bearings for quiet operation. Durable 1.75-inch tubular, epoxy-coated steel frame. This stand has rear cross bar supports.

**TECHNICAL SPECIFICATION:**

Dimensions, W x D (cm / in)	182.9 x 101.6 / 72.0 x 29.0
Maximum Stand Height (cm / in)	85.1 / 34.5
Minimum Stand Height (cm / in)	69.9 / 27.5
Nominal Width (cm / ft)	182.9 / 6.0
Weight (kg / lb)	57.0 / 125.0
Conformance	ADA