



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6769**





### **FLIP-OFF SEAL (YELLOW, 20 mm)**

Plastic w/Aluminum base

**DESCRIPTION:** Seals are a crucial part of a vial system. The selection of an appropriate seal is critical to ensuring container closure integrity, drug protection, user convenience and safety. To deliver consistency and extreme precision, Flip-off seals are manufactured using TruEdge<sup>®</sup> technology production that includes 100% vision inspection.

#### **TECHNICAL SPECIFICATION:**

-  The plastic top flips off easily with pressure from the thumb exposing the aluminum seal with center hole
-  Use on 20 mm OD vial