



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6770





FLIP-OFF SEAL (WILLOW GREEN, 20 mm)

Plastic w/Aluminum base

DESCRIPTION: Seals are a crucial part of a vial system. The selection of an appropriate seal is critical to ensuring container closure integrity, drug protection, user convenience and safety. To deliver consistency and extreme precision, Flip-off seals are manufactured using TruEdge[®] technology production that includes 100% vision inspection.

TECHNICAL SPECIFICATION:

-  The plastic top flips off easily with pressure from the thumb exposing the aluminum seal with center hole
-  Use on 20 mm OD vial