



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 6773





FLIP-OFF SEAL (Red, 13 mm)

Plastic w/Aluminum base

Description: Seals are a crucial part of a vial system. The selection of an appropriate seal is critical to ensuring container closure integrity, drug protection, user convenience and safety. To deliver consistency and extreme precision, Flip-off seals are manufactured using TruEdge[®] technology production that includes 100% vision inspection.

TECHNICAL SPECIFICATION:

-  The plastic top flips off easily with pressure from the thumb exposing the aluminum seal with center hole
-  Use on 13 mm OD vial

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

WWW.MEDISCA.COM