

PRODUCT NO.: 6779



## Mold, Suppository Aluminum 2g (100 Cavity)

**Description:** Aluminum suppository molds produce 2 gm suppositories in a universal shape for rectal or vaginal use. Quick and easy assembly - no tools required. Threaded posts and wingnuts secure the mold in precise alignment.

### Features:

- Dimensions, L x W x H (in): 11.75 x 3 x 1.5
- Weight (lbs): 4

### Maintenance:

Wash mold with soap and water for proper cleaning. Do NOT put in a dishwasher.

USA | Canada | Australia | International

[medisca.com](http://medisca.com)