



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

## PRODUCT NO. 6790











### 2ft TELESCOPING BASE STAND WITH FIXED FEET (LABCONCO)


For XPert / RXPert Filtered Balance Systems

**DESCRIPTION:** Telescoping Base Stands with Fixed Feet may be adjusted to eight height positions prior to installation. Telescoping legs are adjustable in 1" increments from 27.5" to 34.5". Durable 1.75" tubular, epoxy-coated steel frame.

#### FEATURES:

-  Weight: 80.0 lbs
-  Weight metric: 36.0 kg
-  Dimensions: 24.0" w x 29.0" d
-  Dimensions metric: 61.0 x 73.7 cm
-  Conformance: ADA
-  Maximum Stand Height: 34.5" (85.1 cm)
-  Minimum Stand Height: 27.5" (69.9 cm)
-  Nominal Width: 2'

#### RELEVANT PRODUCTS:

-  Applicable Products: Paramount Ductless Enclosures, Protector XVS Ventilation Stations, Purifier Class I Enclosures, Purifier Filtered PCR Enclosures, Purifier HEPA Filtered Enclosures, Purifier Vertical Clean Benches, XPert Balance Enclosures, XPert Filtered Balance Stations, Xpert & RXPert Filtered Balance Systems, XPert Nano Enclosures