



**PRODUCT NO.: 6817**



### **Seals, Aluminum, Flip-Off, 20mm, Blank Top (Blue)**

**Description:** Seals are a crucial part of a vial system. The selection of an appropriate seal is critical to ensuring container closure integrity, drug protection, user convenience and safety. To be used on a 20 mm OD vial.

#### **TECHNICAL SPECIFICATIONS:**

Color	Blue
Inner diameter	20.20 mm $\pm$ 0.20 mm
Thickness	0.20 mm $\pm$ 0.02 mm
Depth	7.30 mm $\pm$ 0.20 mm
Plastic outer diameter	22.20 mm $\pm$ 0.20 mm
Total height	8.80 mm $\pm$ 0.30 mm

USA | CANADA | AUSTRALIA PACIFIC | EUROPE  
[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)