




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO. 6826



CERVICAL SPECIMEN COLLECTION, SWAB W/CERVISOFT, FOAM TIP 7"L, Sterile, Polystyrene handle

DESCRIPTION: This cervical specimen collection device is made with a non-abrasive, soft tip that is gentle on patients and reduces discomfort. The foam head readily releases cervical cells onto slides for testing. Ideal to use for cervical specimen collection, sexual assault kit component, and DNA collection.

FEATURES:

- ✔ 7 inch cytology collection device with non-abrasive cervisoft foam tip and polystyrene shaft
- ✔ Soft, non-abrasive foam tip is gentle on patient
- ✔ Foam head readily releases the endocervical cells onto the microscope slide
- ✔ Sterile and individually wrapped
- ✔ Made in USA

TECHNICAL SPECIFICATION:

- ✔ **Packaging:** 10 boxes of 100 individually wrapped = 1,000/cs
- ✔ **Weight:** 7 lbs.
- ✔ **Overall Length:** 7" (177.8 mm)
- ✔ **Tip Dimension:** 1" Length (25.4mm) x 0.188" Width (4.775mm)
- ✔ **Handle Dimension:** 7.010" Length (178.054mm) x 0.100" Dia. (2.54mm)
- ✔ **Handle:** Polystyrene
- ✔ **Tip Material:** Cervisoft Foam
- ✔ **Sterile:** Sterile