PRODUCT NO. 6833





TAP-N-CLICK MULTI-PURPOSE, SILICONE PAD

Applicator Accessories, 5/pack, Thick silicone material

DESCRIPTION: The Tap-n-CLICK was designed to be a multiple purpose product serving the patient as well as the pharmacy. For the patient, these innovative pads allow for "onehanded" CLICK-Dosing of medication. For the pharmacy, these pads provide a sound dampening surface that assists with the tapping out of air bubbles during filling of medication.

FEATURES:

- M Simple, smart and convenient topical dosing accessories
- Made of thick silicone material
- Easy to tap the Topi-CLICK on the pad to remove air bubbles
- Facilitates one-handed dispensing by using the Pad and turning the Super-CLICKer



For Patients: Click-It



🔚 1.800.932.1039 | 💵 1.800.665.6334 | 📰 1300 786 392 | WWW.MEDISCA.COM