



PRODUCT NO.: 6871



SPILL TRAY, AIR SCIENCES (3ft)

Polypropylene, For Purair Basic Ductless FumeHood (Extra Tall, Reduced Depth), Black

DESCRIPTION: This chemically resistant spill tray is used with the Purair Basic Ductless Fume Hood (6872-01) to contains spills and simplify cleaning procedures.

TECHNICAL SPECIFICATIONS:

- Dimensions, W x D x H (in / cm): 36 x 24 x 1 / 91 x 61 x 2.5
- Weight (lbs / kg): 15 / 6.8

ADDITIONAL INFORMATION:

Applicable Products:

- 6872: 3ft Purair Basic Ductless Fume Hood (Extra Tall, Reduced Depth)