



PRODUCT NO.: 6910



EXHAUST COLLAR, AIR SCIENCES (3ft)
Polypropylene, For PurAir Basic Ductless FumeHood

DESCRIPTION: The exhaust collar allows the enclosure to be thimble-ducted to the outside for additional protection.

TECHNICAL SPECIFICATIONS:

- Dimensions, W x D (in / cm): 32 x 22 / 77.6 X 56.2
- Collar Diameter (in / cm): 8 / 20.3

ADDITIONAL INFORMATION:

Applicable Products:

- 6912-01: 3ft Basic Ductless FumeHood