




**PRODUCT NO.: 6928**



**EXHAUST COLLAR, AIR SCIENCES (4ft)**  
Polypropylene, For PurAir Biological Safety Cabinet

**DESCRIPTION:** The exhaust collar allows the enclosure to be thimble-ducted to the outside for additional protection.

**TECHNICAL SPECIFICATIONS:**

- Dimensions, W x D x H (in / cm): 30.5 x 22.1 x 13.8 / 77.6 X 56.2 X 35.0
- Collar Diameter (in / cm): 8 / 20.3

**ADDITIONAL INFORMATION:**

**Applicable Products:**

- 6930-01: 4ft PurAir Biological Safety Cabinet