











PRODUCT NO.: 6978



Bottle, Amber Glass, 1oz (30mL) w/ Brush Cap

Description: Amber Boston Round Bottles protect contents from UV rays and are ideal for light sensitive products. The black Polypropylene brush caps have a 4 mm diameter PE rod ending in a soft applicator bristled brush. The applicator brushes can be used in applications where the product needs to be brushed on.

FEATURES:

-  Style: Boston Round Bottle
-  Bottle Size, Ø x H (mm): 31 x 79
-  Material: Glass
-  Color: Amber
-  Neck: 20-400
-  Cap: Black Polypropylene Cap w/PE Rod & Soft-Bristled Brush Applicator
-  Rod Size, Ø x H (mm): 4 x 69
-  Volume: 1oz / 30mL

USA | CANADA | AUSTRALIA PACIFIC | EUROPE
WWW.MEDISCA.COM