




**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 6993**



### **OINTMENT TUBE TRAY (U-Heat Seal Solutions)**

Holds 4 x 1 or 2oz Ointment Tubes, HDPE, For 7050-01 Tube Sealer

**DESCRIPTION:** This customized tube tray will hold four 1 or 2 oz plastic ointment tubes while using the 7050-01 Tube Sealer. Please note that this tube tray has been customized to fit the plastic ointment tubes sold by MEDISCA.

#### **ADDITIONAL INFORMATION:**

##### **Applicable Products:**

 7050-01: Ointment Tube Sealer