



PRODUCT NO.: 7045



HEPA FILTER, AIR SCIENCES (4ft) For Purair Basic Ductless FumeHood

DESCRIPTION: Two HEPA filters are required with the purchase of the Purair Basic Ductless Fume Hood (7047-01). The Purair Basic and these HEPA filters will capture all powders and particulates, making the hood to be commonly known as a non-sterile powder containment hood.

TECHNICAL SPECIFICATIONS:

- Dimensions, W x D (in / cm): 13 x 13 / 33 x 33
- Weight (lbs / kg): 5 / 2.3

ADDITIONAL INFORMATION:

Applicable Products:

- 7047-01: 4ft Purair Basic Ductless Fume Hood
- 7046-01: Spillage Tray