



**PRODUCT NO.: 7046**



### **SPILL TRAY, AIR SCIENCES (4ft)**

Polypropylene, For Purair Basic Ductless FumeHood, Black

**DESCRIPTION:** This chemically resistant spill tray is used with the Purair Basic Ductless Fume Hood (7047-01) to contain spills and simplify cleaning procedures.

#### **TECHNICAL SPECIFICATIONS:**

- Dimensions, W x D x H (in / cm): 48 x 27 x 1 / 122 x 68.6 x 2.5
- Weight (lbs / kg): 10 / 4.5

#### **ADDITIONAL INFORMATION:**

##### **Applicable Products:**

- 7047-01: 4ft Purair Basic Ductless Fume Hood
- 7045-01: HEPA filter (2 are required for 1 x 7047-01)