



PRODUCT NO.: 7059



CONVERSION KIT, PROFILLER 1100 (Size 000)

Includes Orienter Base and Locking Plate

DESCRIPTION: A conversion kit is required when switching to a size 000 on the ProFiller 1100. The Conversion Kit includes an orienter base and capsule locking plate. This conversion kit needs to be ordered with size 000 change parts to successfully fill size 000 capsules.

ADDITIONAL INFORMATION:

Related Products:

-  7060: Change Parts, Size 000
-  7058: Tamper, 40 Pin, Size 000