



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7064








MINI STOPPERING CHAMBER (Labconco)

Includes 17" Clear PVC Tubing, Holds up to 48 x 2mL Serum Vials or 8 x 30mL Serum Vials

DESCRIPTION: The Mini Stoppering Chamber attaches to a port on a Labconco Drying Chamber or Manifold (sold separately) to stopper small vials (up to 30ml each). Includes 17" of clear PVC tubing. Can stopper up to 48 serum vials (2ml) or 8 serum vials (30ml) while maintaining vacuum. Spacers, tubing, and adapter for connection to the port are included. Vials and stoppers sold separately.

TECHNICAL SPECIFICATIONS:

-  Weight: 15.0 lbs
-  Weight metric: 6.8 kg
-  Dimensions: 12" w x 12" d x 10" h
-  Dimensions metric: 30.5cm x 30.5cm x 25.4cm
-  Applicable Products: Drying Chambers, Manifolds

ADDITIONAL INFORMATION:

Applicable Products:

-  7063-01: Vacuum Pump, 117L/min
-  7065-01: 4-Port Manifold
-  7066-01: 2.5L Benchtop FreeZone Dry System