



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 7092-01**



**BASE STAND w/CASTERS, AIRCLEAN SYSTEMS**  
4ft, For PowderSafe Balance Enclosure Type A

**DESCRIPTION:** Constructed of 14-gauge epoxy coated steel, AirClean® Systems custom carts are designed exclusively for our ductless hoods, workstations and enclosures. All carts feature a two-position movable base bar for wheelchair or lab chair compatibility. Standard height of 34" meets ADA requirements. An optional polypropylene shelf may be added to each cart allowing storage or transportation of supplies.

**ADDITIONAL INFORMATION:**

**Applicable Products:**

-  8167-01: 4ft AirClean PowderSafe Type A