



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7152-01









PACKER ROUND PLASTIC BOTTLE w/SCREW CAP

100mL, 38-400, PET, Wide Mouth w/White Ribbed PP Cap, Amber

EA

DESCRIPTION: Amber Boston Round Bottles protect contents from UV rays and are ideal for light sensitive products. These general use bottles are perfect for liquids. They feature a round body, rounded shoulders and wide screw neck opening. These environmentally sensitive bottles help eliminate waste and help to insure product integrity for long term storage.

FEATURES:

-  Style: Packer Round Bottle
-  Color: Amber
-  Material: PET
-  Neck: 38-400
-  Attached Cap: White Ribbed PP Cap with SturdeeSeal PE Foam Liner
-  Capacity/Volume: 100mL