



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7154-01









PACKER ROUND PLASTIC BOTTLE w/SCREW CAP

500mL, 45-400, PET, Wide Mouth w/White Ribbed PP Cap, Amber
EA

DESCRIPTION: Amber Boston Round Bottles protect contents from UV rays and are ideal for light sensitive products. These general use bottles are perfect for liquids. They feature a round body, rounded shoulders and wide screw neck opening. These environmentally sensitive bottles help eliminate waste and help to insure product integrity for long term storage.

FEATURES:

-  Style: Packer Round Bottle
-  Color: Amber
-  Material: PET
-  Neck: 45-400
-  Attached Cap: White Ribbed PP Cap with SturdeeSeal PE Foam Liner
-  Capacity/Volume: 500mL