




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7270-01



PRE-FILTERS

3ft, For AirClean Type B (AC730C) PowderSafe Balance Enclosure

PK/12

DESCRIPTION: The AirClean Systems PowderSafe ductless enclosures utilize high efficiency, electrostatically charged pre-filters. Pre-filters are designed to capture larger particulate, aid in minimizing clogging of the HEPA filter, and prevent 'particle rain'.