



PRODUCT NO.: 7277






HEPA FILTER

3ft/5ft, For AirClean Type A/Type B PowderSafe Balance & Bulk Powder Handling Enclosures

DESCRIPTION: Replacement HEPA filter to be used for the 3ft PowderSafe Balance Enclosures and 5ft Bulk Powder Handling Enclosure.

ADDITIONAL INFORMATION:

Applicable Products:

-  7454-01: 3ft PowderSafe Balance Enclosure, Type A
-  7564-01: 3ft PowderSafe Balance Enclosure, Type B
-  6947-01: 5ft Bulk Powder Handling Enclosure

For HEPA filter replacement,

7454-01 will require: 1 x 7277-01

7564-01 will require: 1 x 7277-01

6947-01 will require: 4 x 7277-01