



**PRODUCT NO.: 7278**





**PRE-FILTERS (3ft/5ft)**

For AirClean Type A PowderSafe Balance & Bulk Powder Handling Enclosures

**DESCRIPTION:** Replacement pre-filters to be used for the 3ft Type A PowderSafe Balance Enclosure and 5ft Bulk Powder Handling Enclosure.

**ADDITIONAL INFORMATION:**

Applicable Products:

-  7454-01: 3ft PowderSafe Balance Enclosure, Type A
-  6947-01: 5ft Bulk Powder Handling Enclosure

For pre-filter replacement,

7454-01 will require: 1 pre-filter

6947-01 will require: 4 pre-filters