



**MEDISCA®**  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 7453**



## **Balance Accessories Calibration Weight, 0.01g**

**DESCRIPTION:** Can be used as a reference standard in calibrating other weights and is appropriate for calibrating high-precision analytical balances with a readability as low as 0.1 mg to 0.01 mg.

### **FEATURES:**

- ✔ Tolerance: 0.010mg
- ✔ Style: Flat with one end turned up for easy handling
- ✔ Material: 304 Stainless Steel
- ✔ Density: 7.95 g/cm<sup>3</sup> at 20°C
- ✔ Dimensions (Weight Only): 0.184" x 0.153" x 0.080"
- ✔ Dimensions (Weight in Case): 2" x 2.5"
- ✔ Case: Supplied in high quality polycarbonate case with red insert

**USA | CANADA | AUSTRALIA PACIFIC | EUROPE**

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)