



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 7460-01**










## **BOSTON ROUND GLASS BOTTLE w/SCREW CAP**

2oz (60mL), 20-400, Amber Glass w/Black Ribbed Phenolic Screw Cap

PK/144

**DESCRIPTION:** Amber Boston Square Bottles protect contents from UV rays and are ideal for light sensitive products. These general use bottles are perfect for liquids.

### **FEATURES:**

-  Style: Boston Round Bottle
-  Size: 39 x 94 mm
-  Color: Amber
-  Material: Glass
-  Neck: 20-400
-  Attached Cap: Black Ribbed Phenolic Screw Cap w/PolyCone Liner
-  Capacity/Volume: 2oz (60ml)