

## PRODUCT NO.: 7461



## BOSTON ROUND GLASS BOTTLE w/SCREW CAP (1 oz / 30 mL)

Amber, Black Ribbed Phenolic Screw Cap

**DESCRIPTION:** Amber Boston Round Bottles protect contents from UV rays and are ideal for light sensitive products. The general use bottles are perfect for liquids.

## **FEATURES:**

Style: Boston Round Bottle

Attached Cap: Black Ribbed Phenolic Screw Cap w/PolyCone Liner

## **TECHNICAL SPECIFICATION:**

Size (mm / in)	33.0 × 79.0 / 1.3 × 3.1
Capacity / Volume (mL / oz)	30 / 1
Neck Finish	20 – 400
Color	Amber
Material	Glass