




**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 7478-01**






## **USB INTERFACE KIT w/CABLE (A&D)**

For FX, FZ Series Balances

**DESCRIPTION:** The FX and FZ series balances do not come with a USB port upon purchase. 7478-01 will create a USB interface/port on the balance, allowing for USB connectivity. USB cable included.

### **ADDITIONAL INFORMATION:**

#### **Applicable Products:**

-  7593-01: FX-300i
-  8162-01: FZ-2000i
-  8179-01: FZ-1200i