




MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7562



Hood Accessories, Balance Vibration Isolator, Labconco

DESCRIPTION: Isolator minimizes vibration to sensitive balances. Includes a 13.4" w x 20.4" d x 2.7" thick marble slab, four isolator pads, and a type 304 stainless steel protective top cover. Shipping weight is 70 lbs (31.8 kg).

APPLICABLE PRODUCTS:

 XPert & RXPert Filtered Balance Systems, XPert Filtered Balance Stations