



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7584-01



HOT PLATE SUPPORT ROD

15" L x 0.5" Ø, Hard Aluminum

DESCRIPTION: To be used in conjunction with a hook connector (8448-01) and extension clamp (8449-01) to conveniently hold a thermometer above a beaker or flask. Alternatively, these three products can even be used to hold the beaker or flask itself.

TECHNICAL SPECIFICATIONS:

Length (in / cm)	15 / 38
Diameter (in / cm)	0.5 / 1.3
Material	Hard Aluminum