

PRODUCT NO.: 7589-01



WORK SURFACE PLATE

5ft, Black Epoxy Dished, For Labconco XPert Filtered Balance System

DESCRIPTION: Black chemical-resistant work surface is contoured to fit the dimensions of the enclosure and to contain spills.

FEATURES:

- Weight: 125.0 lbs
- Weight metric: 56.0 kg
- Dimensions: 60.0" w x 29.0" d x 1.0" h Dimensions metric: 152 x 74 x 2.5 cm
- Applicable Products: Paramount Ductless Enclosures, Purifier Class I Enclosures, Purifier Filtered PCR Enclosures, Purifier HEPA Filtered Enclosures, Purifier Vertical Clean Benches, XPert Filtered Balance Stations, Xpert Filtered Balance Systems
- Cupsink cutouts: No cupsink cutout Hood, Enclosure or Glove Box Width: 5'

ADDITIONAL INFORMATION:

Applicable Products: 7592-01