



**MEDISCA**<sup>®</sup>  
YOUR TRUSTED PARTNER IN COMPOUNDING

**PRODUCT NO.: 7709-01**



## **HEPA FILTER (AirClean Systems)**

2ft/4ft, For AirClean Type A PowderSafe Balance Enclosure

**DESCRIPTION:** Replacement HEPA filter to be used for the 2ft or 4ft PowderSafe Balance Enclosures.

### **ADDITIONAL INFORMATION:**

#### **Applicable Products:**

-  7354-01: 2ft PowderSafe Balance Enclosure
-  8167-01: 4ft PowderSafe Balance Enclosure

For HEPA filter replacement,

7354-01 will require: 1 x 7709-01

8167-01 will require: 2 x 7709-01