







**PRODUCT NO.: 7788**



**BOSTON ROUND GLASS BOTTLE w/PLASTIC APPLICATOR ROD (1 oz/30mL)**  
20-400, Amber Glass w/Black Screw Cap

**DESCRIPTION:** This amber round glass bottle with plastic applicator rod will provide controlled non-metered dropping. Applicator has blunt end and is permanently attached to cap.

**FEATURES:**

-  Glass: Amber Round
-  Cap: Black Screw Cap w/Plastic Applicator Rod
-  Neck: 20-400
-  Volume: 1 oz / 30mL