

PRODUCT NO.: 7910



Syringe, BD Luer-Lok™ Tip, 10 mL (Sterile)

Description: Sterile syringes feature an integral BD Luer-Lok™ hub, assuring a secure, threaded connection to needles or IV ports. The clear polycarbonate barrel provides a readily visible fluid level. Syringe only. Individually wrapped.

Features:

- Sterile via Radiation, EtO
- Latex-Free
- Graduated Interval (GI): 0.2mL
- Capacity: 10mL

USA | Canada | Australia | International

medisca.com