



MEDISCA[®]
YOUR TRUSTED PARTNER IN COMPOUNDING

PRODUCT NO.: 7912







LUER LOCK SYRINGE, MONOJECT (3mL)

Sterile, GI: 0.1 mL, Clear

DESCRIPTION: Sterile syringes feature an integral BD Luer-Lok hub, assuring a secure, threaded connection to needles or IV ports. The clear polycarbonate barrel provides a readily visible fluid level. Syringe only. Individually wrapped.

TECHNICAL SPECIFICATIONS:

-  Sterile via Radiation, EtO
-  Latex-Free
-  Graduated Interval (GI): 0.1 mL
-  Capacity: 3mL