



**MEDISCA®**  
YOUR TRUSTED PARTNER IN COMPOUNDING



**PRODUCT NO.:7960**



## **Pump Accessories, Tubing 3/8 x 7/8 x 1/4 (for 8159)**

**DESCRIPTION:** This thick-walled PVC vacuum tubing can withstand a full vacuum (29.9 in or 760mm/Hg at room temperature).


### **FEAUTURES:**

-  Flexible tubing
-  Autoclavable

### **TECHNICAL SPECIFICATIONS:**

Maximum Pressure	27 psi at 125°F (52°C), 54 psi at 73°F (23°C),
Inner Diameter (in)	3/8
Outer Diameter (in)	7/8
Overall Length (m/ft.)	3.05/10
Thickness (cm)	0.64
Color	Clear
Material	PVC

### **APPLICABLE PRODUCTS:**

-  8159: Pump, Vacuum Pressure 60psi (Gast)

USA | CANADA | AUSTRALIA PACIFIC | EUROPE

[WWW.MEDISCA.COM](http://WWW.MEDISCA.COM)