



PRODUCT NO.: 7985



CHANGE PARTS, PROFILLER 3700 (Size 00)

DESCRIPTION: Change parts* for one capsule size include: Caps Tray / Body Sheet Set (A, B, C) / Orienter Tray / Locking Plate Spacer

ADDITIONAL INFORMATION:

*Tamper may be required. Please provide your MEDISCA Sales Representative with the FB serial # (located on right side of the Filler Base) and the capsule sizes currently being filled on your current capsule machine to ensure that all pieces are included for proper capsule filling.

Related Products:

-  7490: Tamper, 300 Pin, Size 00