




PRODUCT NO.: 7986




CHANGE PARTS, PROFILLER 3700 (Size 4)

DESCRIPTION: Change parts* for one capsule size include: Caps Tray / Body Sheet Set (A, B, C) / Orienter Tray / Locking Plate Spacer

ADDITIONAL INFORMATION:

*Tamper may be required. Please provide your MEDISCA Sales Representative with the FB serial # (located on right side of the Filler Base) and the capsule sizes currently being filled on your current capsule machine to ensure that all pieces are included for proper capsule filling.

Related Products:

-  7180: Bottom Plate, Size 4-5 – REQUIRED
-  7199: Tamper, 300 Pin, Size 4