



PRODUCT NO.: 8114










RS232 INTERFACE KIT, A&D For EJ Series Balances

DESCRIPTION: Once installed, this interface kit will allow the balance to be connected to the A&D printer (8284-01) via RS-232C connection.

ADDITIONAL INFORMATION:

Applicable Products:

-  8284: DOT MATRIX PRINTER w/CABLE, A&D (For EJ, FX, FZ & HR Series Balances)
-  8309: EJ-200
-  7579: EJ-410
-  7391: EJ-610
-  8293: EJ-1500
-  8285: EJ-3000
-  8368: EJ-6100